

Stainless Steel Warming Cabinets

ICI Healthcare warming cabinets are designed and manufactured to maintain ideal clinical temperatures for blankets, solution bags and solution bottles to aid in the high quality care of patients.



Standard Construction Details:

- Free standing and counter top designs
- Double wall stainless steel construction with a #4 finish
- Temperature range 90°-160° F (32°-71° C)
- One year warranty
- Fahrenheit or Celsius display
- Heating compartments fully insulated for high efficiency
- One to three adjustable perforated stainless shelves (depending upon model)
- Available with stainless steel double paned or tempered glass doors with extruded aluminum frames
- Field reversible left or right swing hinged doors
- Doors have self closing spring loaded hinges
- Magnetic silicone door gasket
- Cord (8 feet) with hospital grade plug
- Front access serviceability
- Leveling feet
- Stainless steel finished backs

Operating & Safety Controls:

- Illuminated digital power & temperature controls for individual chambers
- Adjustable set point temperature range with set points +/- 5° F
- Temperature lock-in capability
- Recirculating fan in each compartment
- Audible and visual over-temperature indicator
- Automatic shut off for temperatures exceeding 180° F (82° C)

This model available for shipment in 5 days

All warmers are ETL and UL61010-1 listed for US. We offer recessed, in wall models as well as a variety of other standard and custom sizes & styles. Please consult your Dealer for further information and assistance.

Caster Kit



Caster Kit Includes:

- (2) swivel locking casters
- (2) swivel casters
- (2) rear 90 degree bumpers
- (2) front side bumpers
- (12) stainless steel screws

Part Number	WC-CAS
-------------	--------



- Exterior dimensions: 30" W x 74" H x 24" D (76.2cm x 188cm x 61cm)
- Upper chamber dimensions: 24 7/8" W x 15 1/4" H x 23 3/8" D (63.2cm x 38.7cm x 59.4cm)
- Lower chamber dimensions: 24 7/8" W x 36 3/8" H x 23 3/8" D (63.2cm x 92.4cm x 59.4cm)
- Two shelves for lower chamber
- 120 Volt, 1240 watts, 10.8 amps
- Heat gain: 4319 BTU/hr
- **Caster Kit Part Number: WC-CAS**

Part Number	Type	Width
HS930-3024	Solid Doors	30" (76.2cm)
HS931-3024*	Framed Glass Doors	30" (76.2cm)

*This item stocked in limited quantities and can be shipped in 5 days



- Exterior dimensions: 30" W x 74" H x 24" D (76.2cm x 188cm x 61cm)
- Chamber dimensions: 24 7/8" W x 62 1/4" H x 23 3/8" D (63.2cm x 158.1cm x 59.4cm)
- Three shelves included
- 120 Volt, 840 watts, 7.0 amps
- Heat gain: 2867 BTU/hr
- **Caster Kit Part Number: WC-CAS**

Part Number	Type	Width
HS920-3024	Solid Door	30" (76.2cm)
HS921-3024	Framed Glass Door	30" (76.2cm)



Bench Top Style

- Exterior dimensions: 30" W x 26" H x 24" D (76.2cm x 66cm x 61cm)
- Chamber dimensions: 24 7/8" W x 15 1/4" H x 23 3/8" D (63.2cm x 38.7cm x 59.4cm)
- 120 Volt, 445 watts, 3.7 amps
- Heat gain: 1515 BTU/hr

Part Number	Type	Width
HS900-3024	Solid Door	30" (76.2cm)
HS901-3024	Framed Glass Door	30" (76.2cm)



Unless otherwise noted warming cabinets can be shipped in 7 weeks.