

# PROCEDURE CARTS

## ACCESSORIES

- ❑ **WORKTOP & SIDE SHELF \***  
# CRTA-SLSHLF-R  
Worktop: ABS 60x40cm (24x16")  
+ Stainless Steel 40x30cm (16x12")



- ❑ **WORKTOP GALLERY TOP**  
# CRTA-TOPPC-3SGT  
Worktop: ABS 60x40cm (24x16")



- ❑ **WORKING SURFACE**  
# CRTA-TOPPC-ST  
Standard for procedure cart  
Worktop: 60x40cm (24x16")



- ❑ **IV POLE**  
# CRTA-IVPOLE-HA  
Height adjustable IV Pole  
Max height 188cm (74")



- ❑ **COMBINATION LOCK**  
# CRTA-MK39/43ELKIT  
Electronic Lock With four-digit user code



- ❑ **KEY LOCK**  
provided as standard



- ❑ **TILT-BIN HOLDER 2 \*\***  
#CATA-TRL-2  
Tilt-Bin Holder for 2 Rows  
Height above cart: 62.1cm (24.8")



- ❑ **TILT-BIN HOLDER 3 \*\***  
#CATA-TRL-3  
Tilt-Bin Holder for 3 Rows  
Height above cart: 86cm (34.4")



- ❑ **ALUMINUM ACCESSORY RAIL**  
#CRTA-SRAIL-MK/NM  
For mounting various accessories on Carts



- ❑ **HAZARD BIN HOLDER**  
#CRTA-HAZB-HL  
Powder coated steel



- ❑ **SIDE CATHETER BOX**  
#CRTA-CATH-HLD  
10.2x10.2x45cm (4x4x18")  
Steel boxes with flexiglass panel




- ❑ **GLOVE HOLDER**  
#CRTA-GLVBX  
For holding disposable glove boxes 24x16.4x9.2cm (9.6x6.5x3.7")



- ❑ **TRIO GLOVE HOLDER**  
#CRTA-GLOVBX-3  
For holding 3 disposable glove boxes 24x16.4x9.2cm (9.6x6.5x3.7")



- ❑ **SIDE BASKET**  
#CRTA-SBASKET  
Wire Side Basket  
44x16x19cm (17.6x6.4x7.6")



- ❑ **WASTE BASKET**  
#CRTA-WSBIN  
Waste Basket, Large with Holder 26x20x36.5cm (10.4x8x14.6") - 10L



- ❑ **PUSH HANDLE**  
#CRTA-SRAIL-MK/NM  
Aluminum 40cm (16")



- ❑ **FOLDABLE FRAME**  
#CRTA-SFRAME-FL  
For side Shelf/side Tray



Should be ordered with 1 or 2 accessory rails

\* Note: Side shelf with gallery top is optional and should be ordered separately

\*\* Bin should be ordered separately



### MULTIPLE DRAWER FRONTS CONFIGURATIONS

Drawer front configuration can be changed instantly without the use of tools i.e. single set of 10cm (4") drawer can be changed with 2 sets of 5cm (2") drawers .

