

# Ordering Guide

S	TYPE	HEIGHT	WIDTH	STYLE	OPTIONS
<b>Cabinet Type</b>					
B	Base (Inset)				
F	Floor (Tall)				
HB	Hanging Base				
OB	Overlay Base				
MB	Mobile Base				
W	Wall				
ASB	Acid Storage Base				
FSB	Flammable Storage				
<b>Base Cabinet Height</b>					
29"	Sitting Height				
32"	ADA Height				
35"	Standing Height				
<b>Wall Cabinet Height</b>					
24" 30" 36" 48"					
<b>Cabinet Width</b>					
12" 15" 18" 21" 24" 30" 36" 42" 48"					
<b>Cabinet Style</b>					
Reference IOPC catalog for door and drawer configurations					
<b>Option Suffix (for base cabinet)</b>					
-18	Depth(contact IOPC for approval)		-PS	1 Pull Out Shelf	
-LK	Locks		-2PS	2 Pull Out Shelves	
-LTK	Toe Kick On Left Side Of Cabinet		-3PS	3 Pull Out Shelves	
-RTK	Toe Kick On Right Side Of Cabinet		-VK	Vent Kit (only for acid cabinet)	
-SS	Stainless Steel, Type 304		-VD	Vertical Divider	
	Stainless Steel, Type 316		-FB	Finished Backs	
	Bread Board(selected models only)		-EC	End Covers	
<b>Option Suffix (for wall cabinet)</b>					
-16	16" Depth		-TG	Tempered Glass	
-LK	Locks		-LG	Laminate (safety) Glass	
-SS	Stainless Steel, Type 304		-THB	Titan Hanging Bracket	
	Stainless Steel, Type 316				
<b>Option Suffix (for mobile cabinet)</b>					
-18	Depth(contact IOPC for approval)		-T	3/4" THK. Metal Top	
-LK	Locks		-E	3/4" THK. Epoxy Top	
-SS	Stainless Steel, Type 304		-P	3/4" THK. Phenolic Top	
	Stainless Steel, Type 316		-2"	Casters (in lieu of standard 3")	
	Bread Board (selected models only)				
<b>Option Suffix (for tall cabinet)</b>					
-LK	Locks		-LG	Laminate (safety) Glass	
-SS	Stainless Steel, Type 304		-WB	Wardrobe Bracket	
	Stainless Steel, Type 316		-PS	1 Pull Out Shelves	
	Tempered Glass				

# Options Menu

## DOOR & DRAWER FRONTS



- Inset steel  
(Standard Default)



- Flush overlay

## DRAWER SLIDES & CASTERS



- 100 lb. full extension  
(Standard Default)
- 100 lb. full extension  
soft close



- 3" Caster (Default)
- 2" Caster

## HINGES & ROLLERS



- Inset Steel
- Stainless steel 5-knuckle  
(Standard Default)




- Flush Overlay
- Stainless steel 5-knuckle  
(Standard Default)

## MAGNETIC CATCH



- Magnetic catch  
(Standard Default)
- Roller catch

## LOCKS



- Master 5-pin  
(Standard Default)

# Options Menu

## PULLS



- Rectangular brushed aluminum



- Recessed aluminum
- Painted recessed aluminum standard colors available



- Aluminum wire pulls
- Stainless steel wire pulls

## LABEL HOLDERS



- None(Standard Default)
- One each door and drawer
- One per pair of doors and each drawer

## BASE UNIT SHELF



- One piece adjustable (Standard Default)

## GLASS FRAMED DOORS



- 3/16" clear float glass (Standard Default)
- 1/4" laminated safety glass
- 1/4" tempered safety glass

## SECURITY PANELS



- Without security panel (Standard Default)
- With security panel