

CABINETS

Medical cabinets with different ISO dimensions to fit your hospital space. Please contact us for more models and information

Specifications:

Side walls: 19mm melamine-coated chipboard + ABS side panels (interior)
 Plinth: 19mm laminated plywood
 Handle: Anodized aluminum
 Edge: 3mm polypropylene, all around rubber seal
 Tambour door: 16mm, polypropylene roller shutter
 Glass door: 3mm reinforced glass
 Solid Door: HPL, light grey - silver
 Hinges: 270 degree opening
 Color: Polar white
 Locking: key lock or code lock



* Depth changes according to door type

** The cassettes are according to ISO standard and fit Pegasus's different storage systems