

# MEDICAL MODULAR CARTS



provides a unique, diversified storage cart product line. A combination of different enclosures, depths, heights and locking systems enables us to provide the optimum solution for any ward in your hospital or day care clinic.

Four solutions are available:

TAMBOUR DOOR

GLASS DOOR

SOLID DOOR

OPEN CART

Baskets are stored horizontally and can be tilted while pulled out



## OVERVIEW

Fully ISO (60X40cm) compatible

All surfaces are made of polymers and therefore resistant to moisture

ABS top cover, PP bottom floor, PVC exterior walls

The entire inside compartment of the cart is easy to clean

Self-drainage floor for quick and easy cleaning

5" (125mm) medical grade casters, non-marking rubber casters.

2x swiveling + 2x swiveling with brakes

Full 360° wrap-around bumper

Improved dense tray capacity due to the use of high capacity interior side tilted panels. Standard 183cm (72") cart accommodates up to 18 units of 5cm (2") deep trays / 11 pcs of 10cm (4") deep trays / 7 pcs of 20cm (8") deep trays

Extended caster pattern for increased stability

ABS side panels can be equipped with small wheels for smoother and easier movement of heavily loaded trays

Designated perforated rails accommodate the use of telescopic slides for support of heavy items. These can be used in combination with the side panels or separately

E Type



U Type

